


Series P - T-8A cavity, Viton, cartridge, and coil seal kit

Model: 990-208-006



Technical Features

The shelf life for this cartridge seal kit is twenty years from date of purchase. This date will be printed on the kit packaging.

Cartridge Installation Data

	U.S. Units	Metric Units
Cavity	T-8/	4
Cartridge Thread	M16	;
Valve Hex Size	7/8 in.	22,2 mm
Valve Installation Torque	20 - 25 lbf ft	27 - 33 Nm

990-208-006

What models can this kit be used on?

RBAN RBAP

Display Components

	Description	Quantity	View
500-101-116	O-Ring	1	Α
500-102-014	O-Ring	1	В
500-151-300	O-Ring	1	С
530-052-300	Glyd Ring	1	С
515-102-014	Backup Ring	1	В
999-208-006	Literature	1	
375-025	Nut	1	

Instructions/Notes

- Install back-up rings and O-rings as shown in views.
- When installing continuous back-up rings, use the same technique and the same care used to install O-rings.
- After installing the GLYD® Ring, size it to restore original shape and to prevent damage during cartridge installation.
- When installing cartridge be sure to torque to recommended value.
- GLYD® Ring is a registered trademark of W. S. Shamban and Co., Los Angeles, CA, USA